

Mouse Monoclonal anti-Androgen Receptor, Clone GM014

60-0087; 60-0087-7 6 mL; 7 mL predilute Antibody, Ready-To-Use
61-0087; 61-0087-2; 61-0087-5 1 mL; 0.2 mL; 0.5 mL Concentrate Antibody
Isotype IgG1
Concentration See container label

Intended Use

For In Vitro Diagnostic Use.

This product is used to qualitatively detect Androgen Receptor in normal and neoplastic formalin fixed paraffin embedded (FFPE) tissue sections in immunohistochemical detection methodology. Interpretation must be made within the context of the patient's clinical history and other diagnostic test by a qualified pathologist.

Description

Androgen Receptor (AR), a transcription factor, mediates the biological actions of androgens, which are essential for differentiation, development, and maintenance of the male reproductive organs.

Reagent provided

This antibody is diluted in 10 mM Phosphate buffered saline (PBS), pH 7.2 containing 1% bovine serum albumin (BSA) and 0.09% sodium azide (NaN₃) as antimicrobial agent.

Precautions

For professional users.

Proper handling of this product as with any product derived from biological sources according to local and applicable regulations.

Sodium azide is a toxic chemical. The concentration in this product is not classified as hazardous, however, the build-ups of NaN₃ may react with lead and copper plumbing to form highly explosive metal azides. Flush the disposed reagent with large volume of water to prevent azide build-up.

Usage

Dilution

60-0087; 60-0087-7: Ready-To-Use
61-0087; 61-0087-2; 61-0087-5: Dilute 1:50-100 before use when using Acu-Stain™ detection system. Optimum dilution factor may vary depending on the specimen and preparation process and should be determined by each individual investigator.

Staining procedure

Incubate this antibody with tissue section for 30-60 minutes at room temperature. Follow the instructions from the selected detection system.

Positive control tissue

Prostate

Epitope retrieval

Citrate Buffer pH 6

Staining pattern

Nucleus

Storage

Store at 2-8°C.

References

1. Jänne OA, et al. Ann Med 1993; 25:83
2. Loda M, et al. Mod Pathol 1994; 7(3):388-391
3. Bayer-Garner IB, et al. Mod Pathol 2000; 13(2):119-122

Symbols



Catalog No.



Batch No.



In Vitro Diagnostic Use



Temperature Range



Use By

