

## Mouse Monoclonal anti-Melanosome, Clone HMB45

60-0042; 60-0042-7 61-0042; 61-0042-2; 61-0042-5

Isotype

6 mL; 7 mL predilute Antibody, Ready-To-Use 1 mL; 0.2 mL; 0.5 mL Concentrate Antibody

Concentration See container label

**Intended Use** For In Vitro Diagnostic Use.

> This product is intended for laboratory use to qualitatively detect Active Melanosome by light microscopy in normal and neoplastic formalin fixed paraffin embedded tissue sections using immunohistochemical (IHC) methodology. Interpretation of any positive or negative staining shall be supported by a proper control and must be made within the context of the patient's clinical

history and other diagnostic test by a qualified pathologist.

This antibody labels immature melanosomes in normal and pathological tissues. Thus the antibody Description

labels fetal and neonatal melanocyte and junctional nevi. Intradermal nevi and normal resting adult melanocytes are negative. Positive staining does not provide distinction between benign and

malignant melanocytic proliferation.

This antibody may be useful for the identification of melanoma and other tumors showing melanocytic differentiation. Differential identification is aided by the results from a panel of antibodies. In general, tumors of epithelial, lymphoid, glial, and mesenchymal origin are negative.

This antibody is purified immunoglobulin diluted in 10 mM Phosphate buffered saline (PBS), pH Reagent provided

7.2 containing 1% bovine serum albumin (BSA) and 0.09% sodium azide (NaN<sub>3</sub>) as antimicrobial

agent.

**Precautions** For professional users

Proper handling of this product as with any product derived from biological sources according to

local and applicable regulations.

Sodium azide is a toxic chemical. The concentration in this product is not classified as hazardous, however, the build-ups of NaN<sub>3</sub> may react with lead and copper plumbing to form highly explosive metal azides. Flush the disposed reagent with large volume of water to prevent azide build-up.

Usage

Dilution 60-0042; 60-0042-7: Ready-To-Use

> 61-0042; 61-0042-2; 61-0042-5: Dilute 1:50 to 1:100 before use when using Acu-Stain<sup>™</sup> detection system. Optimum dilution factor may vary depending on the specimen and preparation process

and should be determined by each individual investigator.

Incubate this antibody with tissue section for 30-60 minutes at room temperature. Follow the Staining procedure

instructions from the selected detection system.

Positive control tissue Melanoma

**Epitope retrieval** Not Required

Staining pattern Cytoplasmic

Storage Store at 2-8°C.

1. Gown AM, et al. Am J Pathol. 1986 May;123(2):195-203. References

2. Kapur RP, et al. J Histochem Cytochem. 1992 Feb;40(2):207-12.

3. Esclamado RM, et al. Am J Surg. 1986 Oct;152(4):376-85

**Symbols** 

REF

LOT

IVD





Catalog No.

Batch No.

In Vitro Diagnostic Use

Temperature Range

Use By

30042 Rev.03







