

Rabbit Polyclonal anti-CD30 Antibody

60-0188-7 61-0188; 61-0188-2; 61-0188-5 Isotype Concentration	7 mL predilute Antibody, Ready-To-Use 1 mL; 0.2 mL; 0.5 mL Concentrate Antibody N/A See container label
Intended Use	For In Vitro Diagnostic Use. This product is used to qualitatively detect CD30 in normal and neoplastic formalin fixed paraffin embedded tissue sections in immunohistochemical detection methodology. Interpretation must be made within the context of the patient's clinical history and other diagnostic test by a qualified pathologist.
Description	This antibody labels anaplastic large cell lymphoma (ALCL) and Reed-Sternberg cells. It is a useful tool for the identification of ALCL and as a secondary marker for Hodgkin's disease. CD30 antigen is also expressed in other abnormal tissues such as embryonal carcinomas and mesenchymal tumors.
Reagent provided	This antibody is diluted in 10 mM Phosphate buffered saline (PBS), pH 7.2 containing 1% bovine serum albumin (BSA) and 0.09% sodium azide (NaN $_3$) as antimicrobial agent.
Precautions	For professional users. Proper handling of this product as with any product derived from biological sources according to local and applicable regulations. Sodium azide is a toxic chemical. The concentration in this product is not classified as hazardous, however the build-ups of NaN ₃ may react with lead and copper plumbing to form highly explosive
Usage	metal azides. Flush the disposed reagent with large volume of water to prevent azide build-up.
Dilution	60-0188-7: Ready-To-Use
	61-0188; 61-0188-2; 61-0188-5: Dilute 1:100-200 before use when using Power-Stain [™] and Acu-Stain [™] detection systems. Optimum dilution factor may vary depending on the specimen and preparation process and should be determined by each individual investigator.
Staining procedure	Incubate this antibody with tissue section for 30-60 minutes at room temperature. Follow the instructions from the selected detection system.
Positive control tissue	ALCL or Hodgkin's Lymphoma
Epitope retrieval	HIER, Citrate pH 6
Staining pattern	Membrane and Cytoplasm
Storage	Store at 2-8°C.
References	 Schwarting R, et al. Blood. 1989 Oct;74(5):1678-89 de Bruin PC, et al. Leukemia. 1995 Oct;9(10):1620-7 Stein H, et al. Blood. 2000 Dec 1;96(12):3681-95
Symbols	REFLOTIVDImage: Catalog No.Batch No.In Vitro Diagnostic UseTemperature RangeUse By

31720 Rev.00





